NETZSCH
HOSZLEN

Leading Thermal Analysis ■

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2. Run a Baseline Measurement

Baseline measurements for DSC are only necessary for the case of verification, c_p measurements and for testing the system. Usually DSC measurements are sample measurements.

	<ul style="list-style-type: none"> - Open the NETZSCH-Proteus group.
	<ul style="list-style-type: none"> - Select the DSC instrument (e.g. DSC 204 F1).
	<ul style="list-style-type: none"> - Select New in the File menu.

	<p>Measurement Definition Setup</p> <ul style="list-style-type: none"> - Check the instrument setup. - The values of the different properties must be checked and changed (if necessary). - Click Forward.
	<p>Measurement Definition Header</p> <ul style="list-style-type: none"> - Select Measurement Type Correction. - Optionally: Define Laboratory, Project, Operator, Date and Material. - Define Identity and Name. - Select the temperature calibration file*. If enthalpy changes are of interest, select the sensitivity calibration file*. If no enthalpies should be evaluated, "will not be used" can be used. <p>* Suitable for the measurement conditions to be used!</p> <ul style="list-style-type: none"> - Define the gas for Purge 1, Purge 2 and Protective. - Click Forward.

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	<p>Measurement Definition</p> <p>Temperature Program</p> <ul style="list-style-type: none"> - Define the temperature program for the measurement. - Here you can activate the gases and define the gas flow (for instruments with MFCs). Here you can also activate the Intracooler. - Click Forward.
	<ul style="list-style-type: none"> - Define folder and file name and save the measurement.

	<p>Measurement Definition</p> <p>Last Items</p> <ul style="list-style-type: none"> - The previously selected folder and filename is shown in the window (to change folder and filename once again click the Select button). - Click Measure.
	<p>Start the Measurement</p> <ul style="list-style-type: none"> - Push Initial Cond. ON (the gas switches are set to the INITIAL start conditions (carried out automatically during preheating/precooling)). - Click Start. - If necessary, check the Instrument Configuration.

3. Measurement of a sample

a) If a baseline should be subtracted

	<ul style="list-style-type: none"> - After the baseline measurement is finished, select Open in the File menu.
	<ul style="list-style-type: none"> - Select the baseline measurement (carried out in Section 2). - Continue with Open.


Measurement Definition

Fast definition

- Select Measurement Type **Correction+Sample**.
- Set the check mark for **Get these properties from just opened measurement** to define the same settings (e.g. Sample name, Identity ...) as for the correction measurement. To change the further settings (as it often is the case for the Temperature Program), click the corresponding tab (e.g. Temperature Program).

ATTENTION!
Please note that the Temperature Program settings must be identical to the temperature program settings of the correction measurement. The temperature program can only be shortened from behind.

- Enter the **sample mass**.
- To **save the measurement parameters** click on **Select** to define a folder and enter a **Filename**.
- Click **Measure**.

Start the Measurement

- Push **Initial Cond. ON** (the gas switches are set to the INITIAL start conditions (carried out automatically during preheating/precooling)).
- Click **Start**.
- If necessary, check the **Instrument Configuration**.


b) if no baseline should be subtracted

- start with **File – new** (analogous procedure to chapter 2); measurement type: **sample**
- start with **File – open** to call up an already existing file and to change its measurement parameters, if necessary

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
NETZSCH

DSC curve evaluated by means of the autonomous analysis *AutoEvaluation*

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
Search Libraries:

Library
<input type="checkbox"/> Alloys Poster NETZSCH
<input type="checkbox"/> Ceramics Poster NETZSCH
<input type="checkbox"/> Ceramics_Inorganics NETZSCH
<input type="checkbox"/> Elements Poster NETZSCH
<input type="checkbox"/> Metals_Alloys NETZSCH
<input type="checkbox"/> Organics_Food_Pharma NETZSCH
<input checked="" type="checkbox"/> Polymers DSC KIMW
<input type="checkbox"/> Polymers NETZSCH
<input type="checkbox"/> Polymers Poster NETZSCH


<https://www.netzsch-thermal-analysis.com/de/produkte-loesungen/software/proteus>

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